

GENE-Fect™ Transfection Reagent

TLC-001	1.2 ml
TLC-001.5	5 x 1.2 ml
TLC-001.10	10X1.2 ml

This product is for laboratory research only and not for diagnostic use

Description

GENE-Fect Transfection Reagent is a novel nano-particle-based transfection reagent for the efficient transfection experiments of broad range of mammalian cells with low cytotoxicity.

Content

GENE-Fect Transfection reagent	1.2 ml
User manual	1 sheet

Storage Information

Store at 2-8 °C

Notice

Optimization of experimental conditions may be necessary to obtain optimal results.


Protocol

Refer to the table below when you set up transfection

Culture vessel	DNA per well (μg)	Gene-Fect per well (μl)	Volume of dilution buffer for both DNA & Gene-Fect(ul)	Total volume complexes added per well (μl)
96 well	0.06	0.18	6	12
48 well	0.125	0.375	12.5	25
24 well	0.25	0.75	25	50
12 well	0.5	1.5	50	100
6 well/ 35 mm	1	3	100	200
6 cm/ flask 25cm ²	2.5	7.5	250	500
10 cm/ flask 75cm ²	5	15	500	1000
14 cm/ flask 175cm ²	10	30	1000	2000

Use the following procedure to transfect DNA into mammalian cells in a 6-well format. For other formats, see scaling up or down Transfections.

Note : Optimization may be necessary.



1. Adherent Cells : One day before transfection, plate $0.5 \sim 2 \times 10^5$ cells in 2 ml of culture medium (serum & antibiotics) so that cells will be 60~80% confluent at the time of transfection.

Suspension Cells : Just prior to preparing complexes plate $4 \sim 8 \times 10^5$ cells in 2 ml of culture medium.

2. For the transfection (6 well or 35 mm plate), prepare complex as follows

- Dilute DNA with 100 ul diluent (DMEM or Opti-MEM or 150 mM NaCl), then mix well
- Dilute GENE-Fect reagent with the identical diluent used for DNA dilution, then mix well
- Combine the diluted DNA (a.) with diluted GENE-Fect Reagent (b) Mix well and incubate the mixture for 15 ~ 30 min at RT

3. Add the combine mixture (c.) to each well or 35 mm plate containing cells in culture medium.

4. Incubate cells at 37 °C in a CO₂ incubator for 18~48 hr prior to testing for the transgene expression.

(optional : change the culture medium with fresh one 8 hr after transfection)